

BL44XU Macromolecular Assemblies (Osaka University)

This beamline is specially designed to collect high quality X-ray diffraction data of biological macromolecular assemblies, e.g. protein complexes, protein-nucleic acid complexes, and viruses.

Area of research

Crystal structure analysis of biological macromolecular assemblies (e.g. membrane protein complexes, protein complexes, protein-nucleic acid complexes, and viruses)

Keywords

Scientific field

Biological macromolecular assemblies, Protein crystallography, Structural biology

Equipment

High precision goniometer, Imaging plate detector, CCD detector, Cryostream cooler

Source and optics

Main beamline optics is a double-crystal monochromator and a horizontal focusing mirror.

| | |
|--|---------------------|
| Insertion device | In-vacuum undulator |
| Undulator period : λ_u | 32 mm |
| Number of period : N_{period} | 140 |
| Tunable range | 9 ~ 16 keV |

X-rays at sample

| | |
|---------------------|--|
| Energy resolution | $\Delta E/E < 2 \times 10^{-4}$ |
| Photon flux | 10^{12} photons/s |
| Beam size (focused) | 80 μm (H) \times 500 μm (V) (FWHM) |

Experimental stations

Two individual single rotation axis (horizontal and vertical) goniometers are available in the experimental hutch.

Goniometer


| | |
|------------------------------|----------------------------|
| Deviation from rotation axis | +/- 2.5 μm |
| Position reproducibility | +/- 0.001 deg. |
| Sample translation range | +/- 30 mm |
| Sample-to-detector distance | 250 ~ 1000 mm |
| Collimation | Quadrant slits |
| Sample camera | CCD camera (\times 250) |

Cryostream cooler


| | |
|------------------------|--------------------------------|
| Temperature range | 90 ~ 300 K (N_2 gas) |
| | 35 ~ 300 K (He gas) |
| Nozzle inside diameter | 7mm ϕ |

Detectors

DIP-6040 (Bruker-AXS) and PX210 (Oxford Instruments, ANL) can be used. DIP6040 has a CCD detector (Bruker-AXS SMART6500) surrounded by independent six imaging plate detector systems.



Schematic View of Beamline



| | FULL RESOLUTION | 2 × 2 BINNING |
|-------------------------------|--------------------|------------------|
| Pixels | 3072 × 3072 | 1536 × 1536 |
| Active pixels | 3000 × 3000 | 1500 × 1500 |
| Pixel size at input (avg.) | 69 μm | 138 μm |
| Image size | 18 M byte | 4.5 M byte |
| ADC | 16 bit | 16 bit |
| Data transfer time | 5 sec | 5 sec |

DIP6040

Imaging plate detector parts

Number of imaging plates : 6

Size of an imaging plate : 400 mmφ

Readout system : pulse motor control

IP changer : Spiral

ADC : 16 bit

Size of pixel : 100 μm × 100 μm

Dynamic range : > 10⁶

Sensitivity : < 1.0 photon/pixel

Partial reading : Available (< 400 mmφ)

Readout time : 400 mmφ : 120 frame/hr.

CCD part (Smart-6500)

Sensitive area : 165 mmφ

Number of CCDs : 1

Demagnification of fiberoptic tapers : 2.75 : 1

Conversion gain (12 keV) : 18 e⁻/x-ray photons

Pixels : 2048 × 2048

Pixel size at input (avg.) : 80 μm

ADC : 16 bit

Data transfer time : 10 sec.

PX210

Sensitive area : 210 × 210 mm²

Number of CCDs : 3 × 3

Gap between each modules : 280 μm

Demagnification of fiberoptic tapers : 2.86 : 1

Conversion gain (12 keV) : 9 e⁻/x-ray photons

Main data processing, structure determination software

MOSFLM

HKL2000 (DIP6040, SMART6500)

CCP4

Contact information

Atsushi NAKAGAWA

Research Center for Structural and Functional
Proteomics

Institute for Protein Research, Osaka University
3-2 Yamadaoka, Suita, Osaka, 565-0871

Phone: +81-(0)6-6879-4313

Fax : +81-(0)6-6879-4313

e-mail: atsushi@protein.osaka-u.ac.jp

Eiki YAMASHITA

Ring Buildg. D29, SPring-8

Research Center for Structural and Functional
Proteomics

Institute for Protein Research, Osaka University
3-2 Yamadaoka, Suita, Osaka, 565-0871

Phone & Fax: +81-(0)791-58-1814

e-mail : eiki@protein.osaka-u.ac.jp