

## BL47XU R&D (1)

An X-ray undulator beamline is allocated to the first 'R&D beamline' which is to serve both scientific and technical R&Ds for the novel utilization of X-ray undulator beam, especially in the hard X-ray region. The missions are R&Ds for novel optics, various imaging techniques and their applications.

Double crystals are both cooled by liquid nitrogen.

### Area of research

R&D

Microtomography

Microbeam

### Keywords

*Scientific field*

R&D, Microimaging, Microtomography, Microbeam

*Equipment*

Free-spaced experimental hutch, Beam monitor, Fresnel zone plate, Silicon drift detector, Clean booth


### Source and optics

#### X-rays at sample

Energy range	5.3 ~ 37.8 keV
Energy resolution	$\Delta E/E \sim 2 \times 10^{-4}$ (@E = 18 keV)
Photon flux	$4 \times 10^{13}$ ph/s (@E = 18 keV, @100 mA)
Beam size	1.2 mm (H) $\times$ 0.3 mm (V)
(FWHM)	(@E = 18 keV, @40 m from source)

### Experimental stations

Two hutches (experimental hutch 1 and 2) are available. Location of these hutches is shown in Figure 1. A clean booth for the sample preparation is located outside of the experimental hutch 2.






Fig.1. Hutch layout of beamline

### Experimental hutches 1 (EH1)

Inside the EH1 is free-space for various uses. The cubic capacity of the hutch is about 3.5 m × .5 m × 3.3 m (height) and the size of the door is about 1.96 m (width) × 2.28 m (height). Users can install every optical equipment for their purpose. A versatile experimental table (1 m × 1 m × 0.9 m (height)) and a pulse-motor controlled X-Y-Z stages are also available. The resolution of the each stage is 2 μm/pulse. When the EH1 is not used, a helium path is installed in order to prevent the decrease of X-ray flux for the EH2.

### Experimental hutches 2 (EH2)

There are two experimental tables. The schematic diagram of the arrangement of two tables is shown in Fig.1. On the upstream table, there are two X-Y-Z stages. On the downstream table, there are four X-Y-Z stages. These stages are driven by each pulse motor. For example, the resolution of the X-Y-Z stages is 2 μm and 0.08 μm/pulse in the Y-Z direction and X direction respectively at full step mode. For general purpose, we have 4-jaw slit, solid-state detector, scintillation counter, ion chamber, He-Ne laser, X-ray cameras (middle resolution: pixel size 6 μm and high resolution; pixel size 0.5 μm), and data acquisition system. A pinhole chamber (pinholes and slits are switchable in vacuum at a distance of 40 m from light source) has been installed in the optical hutch.

### X-ray microtomography

Upstream table is normally used for the X-ray microtomography system consists of sample stage and image detector stage. Air-bearing type sample rotation stage is used. The wobbling accuracy of the sample rotation stage is less than 0.2 μm. A visible-light converted type image detector (called as Beam Monitor) is used as data acquisition. The spatial resolution is about 1 μm. The maximum field of view

is about 1 mm × 1 mm. It takes about 2-3 hours for one CT scan.

### Microbeam (Scanning microscopy)

Downstream table is normally used for the microbeam experiment. A Fresnel zone plate is used as the focusing device. Probe size is 0.5 μm (vertical) × 0.3 μm (horizontal), and the focal depth is about 1 mm. Useful energy range is 7~15 keV. X-ray flux of the probe is about 10<sup>9</sup> photons/s in maximum. Ionization chambers and silicon drift detector (SDD) are prepared as the detector.




Fig.2. X-ray microtomography setup




Fig.3. X-ray scanning microscopy setup

### Sample data 1

The computer tomographic (CT) images of Sn-Pb eutectic alloys with different growth rates are shown. The micrometer-resolution structures were successfully observed with the system. Bright and dark areas correspond to Pb-rich


and Sn-rich phases, respectively. These figures also show that the system had spatial resolution of at least a few microns. X-ray energy was 29.1 keV. The number of projections was 360 for 180 degree rotation. The exposure time for one projection was 2.0 second. The slice thickness was 0.2  $\mu\text{m}$ .



### Sample data 2

Figures (a) and (b) show intensity profiles of focal probe of X-ray microbeam in vertical direction and horizontal direction, respectively. X-ray energy; 8 keV. Focusing device; Fresnel zone plate (NTT-AT), material; Ta, thickness; 1  $\mu\text{m}$ , outermost zone width; 250 nm.

Figure (c) shows a scanning X-ray microscope image of a Ta test pattern (scan pitch 1  $\mu\text{m}$ ), and its finer scan image (scan pitch 62.5 nm) is shown in (d). Image of 0.2  $\mu\text{m}$  space is clearly resolved.



### Contact information

**Akihisa TAKEUCHI**

SPring-8 / JASRI

1-1-1 Kouto, Mikazuki-cho, Sayo-gun, Hyogo 679-5198

Phone : +81- (0)791-58-0833

Fax : +81- (0)791-58-0830

e-mail : take@spring8.or.jp

**Mitsuhiro AWAJI**

SPring-8 / JASRI

1-1-1 Kouto, Mikazuki-cho, Sayo-gun, Hyogo 679-5198

Phone : +81- (0)791-58-0833

Fax : +81- (0)791-58-0830

e-mail : awaji@spring8.or.jp